whiteClass

DENTSCARE LTDA - Av. Edgar Nelson Meister, 474 - Distrito Industrial - 89219-501 - Joinville - SC - Authorization of Functioning MS P5X44Y0XX28 - CNPJ/Tax ID: 05.106.945/0001-06 - BRAZILIAN INDUSTRY - Registration at ANVISA Nº 80172310039 Technical Authority: Friedrich Georg Mittelstädt - CRQ: 13100147-SC - Brand: FGM®

Rev: 01

Irritante

Irritante

Irritant

6% - 7.5% e 10%

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**At-Home Dental Bleaching System**

**For Professional Use Only**

*Read carefully all the information in this product information sheet before using the product. Keep it for reference until total consumption of the product.*

**Product Description**

White Class with calcium is a dentist-controlled and dentist-supervised at-home bleaching system, based on 4%, 6%, 7½% and 10% hydrogen peroxide gel. It has a neutral pH, high water content, potassium nitrate and sodium fluoride as desensitizers. Compared to the precursor White Class, calcium was added with the aim to minimize enamel demineralization.

The formulation of the White Class with calcium at home bleaching system is based on its precursor, White Class. Therefore, this new

bleaching system possess the same basic properties of its precursor: it is a highly viscous carboxol gel, with a pH level close to neutral (6 to 7), high water content and potassium nitrate and sodium fluoride as desensitizers.

**Presentation**

Contains four syringes with 3g of gel each, 04 tips for application of the gel, 01 dental tray case and instructions for the patient and professional.

Bulk: Fifty syringes with 3g of gel each.

**Basic Composition**

Available in 4%, 6%, 7½% and 10% hydrogen peroxide gel. Neutralized carboxol, potassium nitrate, sodium fluoride, calcium gluconate, stabilizer, deionized water and surfactant.

**Indications**

It is a dentist-supervised product for professional use only for vital teeth bleaching. White Class with calcium was developed for at-home vital tooth bleaching and should be indicated and performed under control and supervision of a dentist. The product is indicated in the following cases:

Physiological darkening of vital teeth due to age.

Very chromatic vital teeth (yellow or darker appearance).

Vital discolored teeth due to incorporation of pigments from some substances such as coffee, tobacco, tea, etc.

Questionable to moderate fluorosis stains.

Low severity of tetracycline stained teeth.

**Precautions and Contra-indications**

The product should be indicated and used only under supervision and control of a dentist.

The product is not recommended for patients who have recently undergone gingival surgery, patients with periodontal disease, sensitive teeth, exposed roots, dental caries, ulcers, poor general and oral health conditions or some allergic reaction to any of the components of the product. These conditions should be treated before the beginning of the bleaching protocol.

Patients with sensitive teeth before the beginning of the bleaching procedure should be treated beforehand.

Teeth with cracks or micro-cracks may present severe teeth sensitivity

or intermittent hypersensitivity, such as stabbing pain. Try to identify and treat these cases previously.

After brushing the teeth, it is recommended to delay the bleaching tray placement with the gel for 30 min. Gingival tissue sensitized by the dental paste or by the mechanical brushing can be worsened by the contact with the hydrogen peroxide.

Dental bleaching trays with cutting edges might cause lesions in the gingival tissue and this can be worsened by the hydrogen peroxide action.

The use during pregnancy or lactation is not recommended.

**Side Effects**

Some patients might experience tooth sensitivity during and after the product application. The intensity of this dental sensitivity is usually mild to moderate and disappears in 1 or 2 days after the treatment interruption. In case the patient experience severe tooth sensitivity, the bleaching treatment should be discontinued and the tooth sensitivity treated before continuing with the treatment. Some patients might have some sore throat due to the swallowing of higher amount of bleaching gel. This occurrence can be prevented by controlling the amount of bleaching gel that the patient place inside the bleaching tray. Non-visible existing hypoplasia or hypocalcification white stains may become more evident during bleaching.

Patients with sensitive or without soft tissue health may experience some gingival sensitivity due to the irritation caused by the contact with hydrogen peroxide.

**Instructions for Use**

Remove the syringe cap and twist the application tip at the bleaching gel syringe.

Place a small drop of the gel in into every compartment of the bleaching tray for the teeth undergoing treatment.

Seat the tray, with the gel, around the teeth with a slight pressure with the aim to distribute gel around all dental surfaces.

Wipe away excess gel in mouth with a dental brush or with the fingers.

The daytime treatment guidelines for the following concentrations are:

4% White Class with calcium: 2 hour/day;

6% White Class with calcium: 1½ hour/day;

7½% White Class with calcium: 1 hour/day;

10½% White Class with calcium: 30 minutes once a day.

Note: Customized treatments can be used according to the etiology of the stain, depending on patient's needs and age, level of sensitivity, etc.

During the product use, follow strictly the dentist's instructions.

Rinse the bleaching trays with water before storing them and before filling them in with the bleaching gel.

Store the bleaching gel in a cold place. Avoid prolonged exposure to light.

**Patient Care**

Some patient's precautions should be adopted before the beginning of the treatment:

1. Examine the patient thoroughly and diagnose the likely cause of the teeth discoloration. Observe in the clinical exam the presence of failed restorations, areas of dental caries, periodontal disease and/or other related problems that might induce tooth sensitivity or pain during the vital bleaching protocol. A periapical radiograph, usually taken under clinical examinations, might be useful before the beginning of any restorative/esthetic procedure.

2. Inform patients about the level of whitening that he is likely to experience during the bleaching treatment. Explain that existing crowns, tooth colored fillings, etc., do not whiten and they will have to be replaced approximately 7 days after the end of the bleaching therapy. Also it is important to report that non-vital teeth do not respond as well as vital ones and some specific bleaching treatment might be required (Whiteness Super Endo or Whiteness Perborato - FGM).

3. Tetracycline stained teeth might not bleach in the same intensity than vital normal teeth. Grayish teeth, even though those not caused by tetracycline use, may not respond well to the bleaching procedure.

4. A slight discoloration may appear in the tray following bleaching adjacent to amalgam fillings (due to superficial oxidation).

5. It is recommended to take pictures before and after the bleaching procedure to observe the results. It is also recommended to bleach one arch at a time in a way that the other arch will serve as reference.

6. Show patients how to load the bleaching tray. Wipe away excess gel in the mouth with a dental brush or with the patient's finger.

7. The prolonged contact of the bleaching gel with the gingival tissue may result in gingival irritation. Deliver to your patient the patient's instruction sheet and emphasize the importance of reading and

understanding it.

8. Inform patients that they might experience some sensitivity in case their teeth are exposed to temperature variations, mainly to cold. In case the tooth discomfort persists, employ the 0.2% Desensibilize KF - FGM (mild to moderate sensitivity) or 2% Desensibilize KF - FGM (in office application from moderate to severe sensitivity) which are desensitizers agents developed exclusively for tooth bleaching.

**Note:** In case the sensitivity is due to exposure of the root by gingival recession or other adverse factors, protect the area or cut out the bleaching tray in order to avoid a prolonged contact of the gel with the sensitive area. However, one should expect that these areas will not bleach with the same intensity than the other regions.

**Selecting the Concentration**

The White Class with Calcium concentrations lead to the same bleaching level. Some factors should be observed to choose one of the concentrations correctly, including the etiology and the level of dental pigmentation, the level of dental hypersensitivity of the patient, the amount of time available for daily use of the tray, etc.

**Conservation and Storage**

Store the product at temperature between 5°C and 25°C / 41°F and 77°F. Do not freeze the product. Keep out of direct sunlight.

**Warnings**

Do not use the product after the expiration date. Discard the product following the legislation of your country. Do not reuse the empty packing. Keep the product out of children's reach.

*This material has been developed solely for use in dentistry. Processing should be carried out strictly according to the instructions for Use. Liability cannot be accepted for damages resulting from failure to observe the Instructions or the stipulated area of application. Besides that, the user is responsible for testing the material for its suitability and use for any purpose not explicitly stated in the Instructions. Descriptions and data constitute no warranty of attributes and are not binding.*