maxxion^C

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Glass Ionomer Cement

Professional Use Only

Read all the information of these directions carefully before using the product. Keep them for consultation, at least until the product has been totally used or there is no more interaction of the product with your last patient.

Description

Maxxion C is a self-curing glass ionomer cement used for luting purposes. The product consists of polycarboxylic acid and glass ionomer, adequately micronized in order to produce a paste with excellent viscosity and the capability to form homogeneous thin films. The main characteristics of this material consist of its fluidity (better settlement of prosthetic components and orthodontic bands), elevated bond strength and the capability of fluoride release. The cement is also radiopaque, which aids on radiographic control.

Forms of presentation

Package containing: 1 bottle containing 15g of powdered cement, 1 bottle containing 10g of liquid, 1 dose dispenser for the powder and 1 mixing pad.

Individual liquid package: 1 bottle of liquid with 10g.

Individual powder package: 1 bottle of cement powder with 15g and 1 powder dispenser.

Basic Composition

Powder: Glass of Fluoro aluminium silicate and Strontium.

Liquid: Polyacrylic Acid and Distilled Water.

Product's Indications

It is indicated for permanent luting of:

crowns, bridges, inlays, onlays (for better adhesive strength, use resin cement. Ie: Allcem , Allcem Core)

Orthodontic bands and other devices (Ie: Haas expander).

For cavity lining under dental restorations.

Note: The cement is compatible with different surfaces: all alloys and acrylic resins.

Precautions and Contraindications

Only for Professional use.

Do not apply the ionomer cement on pulpar exposure. In deep cavities, consider the application of a proper material before the cementation (ie: calcium hydroxide). The adequate proportion of the phases (powder and liquid) must be respected, in a way that the excess of liquid may cause pulpar damage.

Avoid the use of Maxxion C in patients with historical of allergy to glass ionomer cements. In case of allergy, suspend the use immediately and advice the patient to search for medical support. Do not allow the contact of the product with soft tissues or skin. In case of contact, immediately rinse the area with water in abundance.

Avoid contact of the product with the eyes. In case of contact, immediately rinse the area and seek medical support.

Do not mixture the powder or liquid that compose this product with any other glass ionomer product.

Avoid the contact of this product with others that contains Eugenol, since this kind of interaction may affect the cure of glass ionomers.

Side Effects

The product's application directly on pulpar exposure may cause irritation or even irreversible damage (pulp necrosis);

The contact of Maxxion C with products such as eugenol, phenol and P-monochlorophenol, tends to reduce the adhesive strength, therefore it may not be used in this kind of association.

Instructions for Use

Before initiating the treatment read carefully the Warning, Precautions, Contra-Indications and possible Side Effects.

* **Before the cementation process:** Perform isolation of the operatory field:

1. The operatory field should be maintained free of contamination by saliva and / or blood.

Use rubber dam or cotton rolls according to the case.

2. Clean carefully all surfaces: enamel, dentin and the prosthesis, and also dry it correctly.

Avoid excessive drying of the teeth because it could cause post-operative pain to the patient.

Be sure that there aren't any residues of other materials that may damage the efficacy of Maxxion C.

3. For prosthetic pieces it is recommended the use of aluminum oxide blast on the internal surface. This procedure elevates the mechanical retention of the prosthesis. After blasted, the prosthesis should be cleaned before cementation.

Protection of dentin-pulp complex

Avoid applying the glass ionomer cement on the pulp exposure or in too deep cavities.

Porportion and Mixture

1. The recommended proportion for mixture is of 1 level measure of powder to 1 drop of liquid.

2. Shake the bottle before collecting the powder dose. Take a dose of powdered glass ionomer and place it on the mixing pad.

3. Dispense the drop of liquid besides (separated) the powder with the flask on vertical position. Be sure that there aren't bubbles on the drop.

4. Close the flasks before the mixture to prevent the contents to absorb air humidity.

5. The mixture should be conducted on a clean mixing pad or glass plate, with plastic or metallic spatula. Initially, mix half of the powder content with the liquid and incorporate the rest gradually until obtaining a homogeneous mass.

6. Verify the consistency and homogeneity of the mixture and follow the product's application. Avoid exceeding the estimated mixture time, applying the material while it's bright, in order to obtain adequate adhesiveness. Next, there are estimated times on the use process of the product:

Mixing time (agglutination): 30 seconds.

Work time (including agglutination time): 2 minutes.

Set time (from the begginig of mixture): 3-4 minutes.

Note: The estimated times may vary according to work temperature (mouth and ambient temperature).

7.Apply a thin layer of cement in the prosthesis to be cemented. Put the prosthesis in position and make moderate pressure to promote adequate settlement and flow of the cement.

8.Remove eventual cement excess with exploratory probe when it is still on plastic phase. This procedure could be repeated from 5 to 7 minutes after the initial mixture, if necessary.

Preservation and Storage

Maxxion C should be stored in temperature between 15°-30°C/59°-86°F.

Protect the powder against humidity and do not store it in the refrigerator.

Do not store the liquid in refrigerator.

Keep the product bottles in their original package always well sealed. Do not use the product after the expiry date which is 2 years from the date of manufacture..

Warning

Do not use the product after expiry date. In order to discard the product, follow the legislation of your country. Keep out of reach of children.

This material was only manufactured for dental use and must be handled according to the instructions for use. The manufacturer is not responsible for damage caused by other uses or incorrect handling. Besides this, the user is obliged to prove, before use and under his or her responsibility that this material is compatible with the use desired, mainly when this use is not indicated in the instructions for use. Data descriptions do not constitute any type of guarantee and so are not in way linked.

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