Opallis
FLOW

Rev.: 05

BR | EN | ES

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EN Instructions for use

Flowable light-curing composite resin.

Professional use only

Read all the information in this manual carefully before using the product. Keep the manual for later reference at least until the product is finished or until there is no more interaction between the product and its last patient.

Product description

Opallis Flow is a light-curing microhybrid composite resin of medium viscosity used for the restoration of slightly invasive preparations, pit and fissure sealant, base/lining underneath direct restorations, in tunnel-type preparations, radiopaque lining of cavities, repairs of enamel defects, restorations of primary teeth, repairs in composite resin, bonding of tooth fragments, class I, III and V restorations, non-carious cervical lesions. Opallis Flow can be used separately or with Opallis and most of the light-curing composites.

Opallis Flow is light-cured with blue wavelength light between 400 and 500nm (blue light) and with at least 450 mW/cm².

Presentation form

1 syringe containing 2g and 5 application tips, available in the shades: T (Transparent), A1, A2, A3, OA3.5 (Opaque), B0,5, A0,5, OP (Extra-Opaque)

Basic composition

Methacrylic monomers such as Bis (GMA) monomer, Bis (EMA) e TEGDMA, canforquinone, co-initiators, conservants, pigments and 72 % of silanized inorganic load composed of microparticles of barium-aluminum silicate and nanoparticled silicone dioxide with the particles in sizes between 0,05 and 5,0 microns.

Product indication

· Restoration of slightly invasive preparations.

· Pit and fissure sealant.

· Base/lining underneath direct restorations.

· Class I, III and V restorations.

· Tunnel-type preparation.

· Repair of enamel defects.

· Restorations of primary teeth.

· Bonding of tooth fragments.

· Repairs in composite resin.

· Non-carious cervical lesions.

· Planning of preparation walls.

· Restorations of incipient caries.

Precautions and Contraindications

· Use Opallis Flow according to these instructions only.

· The product is indicated for professional use only.

· Avoid using Opallis Flow with patients with a history of allergy to composite acrylates-based resins. In case of allergy, interrupt and ask the patient to seek medical assistance.

· To eliminate the risk of irritation, the areas close to the pulp should receive the appropriate protection (e.g. application of calcium hydroxide).

· To avoid allergic reactions to the methacrylate monomer, the constant use of protective gloves is recommended. If Opallis Flow gets in contact with the hands, wash them abundantly with water.

· Avoid contact with the eyes. In the case of contact with the eyes, wash them abundantly with running water and contact a doctor immediately.

· Use protective goggles when utilizing a light-curing device.

· Opallis Flow is contraindicated when it is not possible to establish a dry operative area or to follow the described technique correctly.

· Opallis Flow should not be ingested or inhaled. If ingested or inhaled, seek medical assistance immediately.

· Keep Opallis Flow out of patients and children's reach.

· Place the syringe lid back after each use.

· Do not reuse the application tip.

Side Effects

The product can cause some irritation level if in contact with the dental pulp and / or very close to it. In cases of allergic reactions to the product, interrupt its use.

Opallis Flow is contraindicated for patients allergic or hypersensitive to methacrylate and related monomers and it should not be used in such cases.

Instructions

- Application in direct restorations:

1. Shade taking: take the shade with the teeth well hydrated.

2. Preparation of the cavity: the preparation of the cavity should be accomplished according to the norms of the adhesive technique, removing only the carious tissue. there is no need for additional retentions.

3. Pulp protection: cover only the areas very close to the pulp (e.g. calcium hydroxide).

4. Etching of enamel and/or dentin: apply phosphoric acid gel Condac 37 (FGM) on the enamel and/or dentin and let it act for 15 seconds. Soon after, wash abundantly with water and dry with oil free air jets. Avoid the dehydration of the dentin.

5. Adhesive agent: apply the adhesive system Ambar (FGM) according to the manufacturer's instructions. If using a self-etching adhesive Ambar Universal (FGM) skip step 4.

6. Application of Opallis Flow: apply the material directly in the cavity with the application tip. Each applied layer should not be thicker than 2 mm. Light-cure each layer for 40 seconds.

7. Finish and polishing: Diamond or Diamond Flex felt disks with the aid of Diamond polishing pastes.

· For pit and fissure sealant:

1. Make the total isolation and the acid attack of the pit and fissure surface, wash and dry.

2. Apply Opallis Flow with the aid of an application tip. Do not allow the product to spread to regions of occlusal contact.

3. Light-cure for 40 seconds. Accomplish the finishing and the

polishing recommend for light-cure resins.

Conservation and storage

Keep the product in a cool place with the packaging always well closed and protected. Protect form the incidence of direct sunlight. Exposure to light causes premature light-curing. Store the product in temperatures from 5°C to 30°C/ 41°F to 86°F. Do not freeze the product. Do not expose the product to high temperatures.

Expiration date

Printed on the product packaging.

Warnings

Do not utilize the product after the expiry date. When disposing of the product, follow the legislation of your country. Keep out of children's reach. The contact of non-light cured Opallis Flow with the skin, mucous membranes, and eyes should be avoided. Before light-curing, the product can cause a slightly irritating effect and sensitivity to methacrylates.

This material was manufactured for dental use only and must be manipulated according to the manufacturer's instructions of use. The manufacturer is not responsible for any damage caused by misuse or improper manipulation of the product. The user is obligated to prove (before using the product and under his/her responsibility) that this product is compatible with the desired utilization, especially when such utilization is not indicated in the instructions of use. Data description does not constitute any type of warranty and therefore is not binding.

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