DENTSCARE LTDA - Av. Edgar Nelson Meister, 474 - Distrito Industrial - 89219-501 - Joinville - SC - Brasil / Brazil. Autorização de funcionamento / Authorization of functioning: MS P5X44XYXX28 - CNPJ / Tax ID: 05.106.945/0001-06 INDÚSTRIA BRASILEIRA / BRAZILIAN INDUSTRY - Responsável Técnico / Technical Authority: Friedrich Georg Mittelstädt CRQ: 13100147 - SC - Registro na ANVISA Registration at ANVISA nº 80172310079 - Marca / Brand: FGM® / Fabricado no Brasil / Made in Brazil

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Composito polimerizado para CAD/CAM

Somente uso profissional

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EN Instructions for use

Cured composite for CAD/CAM

For professional use only
Read carefully all the information in this manual. Save it for further consultation, at least, until the product lasts and/or until there is no more interaction between the product and its last patient.

Product description

Brava Block is a nanohybrid composite, cured at ideal conditions in the form of blocks for milling in automated lathes. The composite presents particles ranging from 40nm to 5.0µm, total inorganic load content in weight of 72% to 82% (52% to 60% in volume) in a highly cured resin matrix. The composite block has a conversion degree of approximately 86%, significantly higher than the majority of the conventional composites (55% to 60%). This feature, along with the adequate distribution of particle sizes, results in a composite of high mechanical properties, mainly its hardness and resistance to wear, which contributes to the application in anterior and posterior teeth. Brava Block is able to provide excellent finishing with smooth and highly shiny surfaces.

Basic composition

Active ingredients: metacrylic monomers, initiator, co-initiator, stabilizers and silane.

Inactive ingredients: vitroceraamic load particles, silica and pigments.

Presentation

Kit: 5 units, 14L.
Low Translucency: A1-LT, A2-LT, A3-LT, A3.5-LT, B1-LT, C2-LT, Bleach.

High Translucency: A1-HT, A2-HT, A3-HT, B1-HT.

Product indication

Brava Block is suitable for temporary or permanent indirect restorations.

- Total unitary crown on tooth
- Inlays and onlays
- Veneers

Precautions and contraindications

For dentistry use only.

Read carefully the Instructions for Use before using the product.

The steps for manufacturing the CAD/CAM systems are described in the use instructions contained in each equipment's manual. The recommendations from manufacturers should be followed;

The milling machine should be set to wet operation while processing the CAD/CAM material. Dry operation is not suitable for Brava Block;

Make sure to leave an internal room for the cement; The piece should be sand blasted in its internal surface, preferably. In case it is not possible, the internal surface can be acid etched with the fluoridric acid at 10% for 60 seconds;

The piece should not be laser-etched;

Do not use hydrogen peroxide solution, EDTA, bicarbonates or astringent products on the surface of the piece;

Do not use temporary cements that contain eugenol;

The piece should not be heated.

Side effects

If the product is used according to the instructions and before the expiry date, no side effect is expected.

Instructions for use

Mechanical preparation of the tooth:

For the piece to have an adequate mechanical performance, it is necessary to follow certain orientations during the dental preparation, according to the procedure to be done:

Inlay and onlay:

The internal axial walls of the preparation should present expansive angles of approximately 5° to allow correct seating of the pieces. Acute corners should be rounded. The occlusal reduction (cusps) should be of 1.5mm and the buccolingual distance

Tabla 1: Inlay



Tabla 2: Onlay


should be at least 1.5mm. (Table 1 e Table 2).

Veneer:

The depth of the preparation on the facial surface varies in the different thirds of the teeth: in the cervical area, it is suggested to be 0.4mm, in the medium third the suggestion is 0.6mm. The incisal reduction recommended is from 0.5mm to 1.5mm. When "enveloping" the incisal border, all angles should be

Tabla 3: Veneer


rounded (Tabla 3).

Total unitary crown:

The axial walls of the preparation must have approximately 5o of inclination and the wear depth must be of 1.0mm in the cervical area and 1.5mm in the medium third. The incisal/occlusal region should be reduced in 1.5mm. All of the preparation should

Tabla 4: Unitary crown

present round angles (Tabla 4).

Occlusal adjustments in the piece and polishing:

Once ready, the piece should be fitted on the tooth before cementation. If any occlusal adjustments are needed, it can be done with fine grit burs. After those adjustments, polishing can be done with rubber or abrasive disks of medium and fine grit, felt disks Diamond (FGM) and the polishing paste Diamond Excel (FGM).

Adhesive preparation of the piece:

1. Sandblast the internal part of the piece using aluminum oxide of up to 50 micrometers, with pressure of around 3 bar until the surface is matted.

As an alternative, the professional can etch the surface with hydrofluoric acid at 10% for 60 seconds.

2. Proceed to the washing of the pieces in an ultrasonic cleaner for 180 seconds.

3. Apply 70% alcohol and dry the surface for 30 seconds.

4. Apply silane (preferably containing MDP) and let it react for 60 seconds. Remove excess with air blasts.

5. Apply the Ambar Universal (FGM) adhesive, which contains MDP for 30 seconds and then apply air blast for 10s to volatilize the solvent and control the thickness of the layer. Then, light cure the whole surface for 20 seconds.

Adhesive preparation of the tooth:

- Using the two-step adhesive Ambar (FGM):

1. Remove any dirt from the tooth preparation. If necessary, use Robinson brush and prophylactic paste. Isolate the tooth, preferably with total dam isolation;

2. Apply the 37% phosphoric acid etchant Condac 37 (FGM) on enamel and dentin for 15 seconds. Then, wash it gently leaving the dentin humid;

3. Apply the adhesive with friction all over the surface of the tooth, for 10 seconds, in two layers. Then, apply air gently to spread the adhesive better and reduce the adhesive layer thickness. Light cure for 20 seconds.

- Using the single-step adhesive Ambar Universal (FGM), which is compatible with FGM's line of resin cements:

1. Remove any dirt from the tooth preparation. If necessary, use Robinson brush and prophylactic paste. Isolate the tooth, preferably with total dam isolation;

2. Apply the adhesive with friction all over the surface of the tooth, for 10 seconds, in two layers. Then, apply air gently to spread the adhesive better and reduce the adhesive layer thickness. Light cure for 20 seconds.

Cementing the piece to the tooth:

1. With the prosthetic piece and the tooth prepared, in an isolated operator field free from fluids, the piece is cemented. For that, manipulate the cement or use the dual-body syringe of Allcem or Allcem Core (FGM) according to the instructions for use of the chosen product;

2. Apply the cement inside the prosthetic piece and place it on the tooth. The excess cement that may overflow must be removed. Light cure each surface of the tooth (labial, palatal/lingual and occlusal) for 20 seconds each, totalizing 60 seconds of light application. After that, keep the field isolated for additional 10 minutes for the cement to cure chemically.

3. Check for proximal excess with dental floss and, if necessary, remove remaining resin/adhesive excess.

Repairing the piece after cementation:

If, at any time, the professional needs to repair a prosthetic piece due to fractures or chipping, the following is suggested:

1. Remove any dirt from the surface of the restoration area. If necessary, use Robinson brush and prophylactic paste. Isolate the tooth, preferably with total dam isolation.

2. Roughen the area to be repaired with fine grit diamond burr, wash and dry the surface;

3. Apply Ambar Universal (FGM) or Prosil (FGM) on the piece, according to these products use instructions;

4. Do the necessary repair by using conventional composite Opallis or Lis (FGM) and light cure the increments according to the product's instructions

for use;

5. Proceed with the finishing and polishing according to the steps of a conventional restoration.

Conservation and storage

Protect from direct sunlight. Store the product between the temperatures of 5°C and 35°C / 41°F and 95°F.

Warnings

- Make sure to use water during the milling process (wet operation). Dry operation is not suitable for Brava Block.

- The product is for professional use only.

- Do not use the product after its expiration date.

- Follow your country's legislation when disposing it.

- Keep out of the reach of children.

- Avoid the use in patients with a history of allergy to the composite resins containing methacrylates. In case of allergy, suspend the use and send the patient to a doctor.

This material was manufactured for dental use only and must be manipulated according to the instructions for use. The manufacturer is not responsible for damages caused for other use or incorrect manipulation. Additionally, the user is obligated to prove, before the use and under his/her responsibility, if this material is compatible with the desirable utilization, mainly in cases in which such utilization is not indicated in this manual. Data description do not constitute any type of warranty and, for this reason, are not binding.

Legenda de símbolos Meaning of symbols Leyenda de los símbolos	
	Data de fabricação Date of manufacture Fecha de fabricación
	Número do lote Batch Number Número del lote
	Prazo de validade Shelf life Fecha de validez
	Limites de temperatura Temperature limits Limites de temperatura
	Não utilizar se a embalagem estiver danificada Do not use if packaging is damaged No lo utilice si el empaque está dañado
	Conservar seco Keep it dry Conservelo seco
	Manter afastado da luz solar Keep protected from sunlight Manténgalo lejos de la luz solar
	Fabricante Manufacturer Fabricante
	Consulte as instruções de utilização Refer to instructions for use Consulte las instrucciones de utilización