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ES Instructions for use

Commercial availability:

The product is commercialized in 1L or 2L containers.

Recommendations:

The OX-19 antiseptic solution, when diluted, is recommended for oral prevention and reduction of the load of microorganisms (including enveloped viruses, such as the one that causes COVID19), present in the aerosol formed during dental treatment procedures. Also, it can be used as an oral colutury helping oral asepsis before and/or after dental procedures, reducing the load of microorganisms present in the oral cavity.

Instructions for use:

- Caution! Concentrated product, must be diluted in water before use in the proportion of 1:20.
- 1- In order to obtain the correct amount for dilution, remove the lid of the dosing chamber, identified with the graded scale;
- 2- Press the bottle in its central part for the liquid to reach the desired dosage for dilution;
- 3- Transfer to the measuring cup the proportion of the concentrated product to the volume desired for the solution. For example, for each liter of water, add 50 mL of concentrated product;
- 4- Transfer to a plastic recipient and add water to the desired final volume;
- 5- Shake for homogenization and it will be ready for use.

Note: the dilution process can be done inside the reservoir (PET bottle) used for the irrigation of the dental equipment. Make sure the recipient contains the final volume corresponding to the correct dilution of the product.

Basic composition

Hydrogen peroxide, citric acid, saccharine, flavoring and purified water.

Precautions

- THE PRODUCT CONTAINS HYDROGEN PEROXIDE.
- USE PERSONAL PROTECTIVE EQUIPMENT (PPE) WHILE HANDLING THE PRODUCT.
- DANGER: CAUSES SEVERE BURNS.
- CAREFUL: IRRITATING FOR THE EYES AND SKIN. DO NOT MIX WITH OTHER PRODUCTS.
- PRECAUTIONS AND CARE IN CASE OF ACCIDENT: do not inhale vapors/aerosols
- IN CASE OF INGESTION: wash mouth. Do not induce vomiting.
- IN CASE OF CONTACT WITH SKIN (OR HAIR) BEFORE DILUTION: remove all contaminated clothes immediately. Wash skin with water/shower. Wash contaminated clothes before wearing them again.
- IN CASE OF CONTACT WITH THE EYES BEFORE DILUTION: carefully wash eyes with water for several minutes. Continue rinsing and seek medical assistance.
- IN CASE OF INHALATION: take the person to a well-ventilated area and keep them in a position that does not make breathing more difficult.
- AFTER DILUTION, IT CONTAINS 0.75% AND DOES NOT REPRESENT ANY RISK TO HEALTH.
- Whenever possible, in case of accident, take the label or the container to the doctor with the patient.
- Do not mix with water in the original container.
- Never reuse the empty container for other purposes.

Contraindications

- This product is not recommended for dental whitening.
- This product is contraindicated for patients with hypersensitivity or allergic to the components of the formula.

Warnings

- Do not dilute or store the solution in metal containers. The contact with metal destabilizes the peroxide and the solution loses its effectiveness in a short time.
- Avoid storing and/or using the diluted solutions for a period longer than 7 days. Try to dilute the product at the time of use, because freshly prepared solutions are in their best condition for use.
- The quality of the water used for the dilution may affect the stability of the solution. Always use filtered, purified or distilled water. The higher the purity of the water, the greater the stability of the peroxide.
- Foam may occur during use.
- Do not ingest the product intentionally. Keep the oral suction device turned on while using the product.
- There may be light tissue irritation in previously sensitized areas, especially in the throat region.

Storage requirements:

- Store product at temperatures between 5° and 25°C.
- KEEP OUT OF REACH OF CHILDREN AND PETS.
- This product must be stored in its original package and closed.
- Keep the product in a dry place away from sunlight.

PROFISSIONAL | PROFESSIONAL | PROFESIONAL

Whiteness do Brasil Indústria Ltda. Rua Dona Francisca - 8300 - Bloco D - Módulo 1 - Joinville - SC - Brasil - CNPJ / Tax ID: 32.256.235/0001-35 **INDÚSTRIA BRASILEIRA / BRAZILIAN INDUSTRY** - Responsável Técnico / Technical Authority: Friedrich Georg Mittelstädt - CRQ 13100147-SC - ANVISA nº: 81867260015 - Marca / Brand: FGM®

Service ao profissional (Brasil) 0800 644 6100
Technical support (Europe) / Servicio al profesional (Americas) +55 47 3441 6100
www.fgm.ind.br | fgm@fgm.ind.br

Legenda de símbolos Meaning of symbols Leyenda de los símbolos	
	Data de fabricação Date of manufacture Fecha de fabricación
	Número do lote Batch Number Número del lote
	Prazo de validade Shelf life Fecha de validez
	Limites de temperatura Temperature limits Límites de temperatura
	Não utilizar se a embalagem estiver danificada Do not use if packaging is damaged No lo utilice si el empaque está dañado
	Conservar seco Keep it dry Conservévalo seco
	Mantenê afastado da luz solar Keep protected from sunlight Manténgalo lejos de la luz solar
	Fabricante Manufacturer Fabricante
	Irritante Irritant Irritating